

Comprehensive solutions for benchtop lab equipment

VWR® PROFESSIONAL HOTPLATE STIRRERS

- Excellent temperature uniformity and powerful stirring
- Large LCD screen presents temperature, speed and time
- Cool-touch, chemical-resistant housing
- Best-in-class five-year warranty

VWR® Professional hotplate stirrers deliver accurate and repeatable results. Temperature range from ambient +5 °C to 500 °C; speed range from 60 to 1600 rpm. Control panel features easy-to-use, intuitive-control knobs to navigate the settings and menu options. Rear housing features an off-centered integral support rod holder with locking knob that works with the optional support rod kit.

| Description | Top plate material | Top plate dimensions | Electrical | Overall dimensions | Weight | Cat. No. |
|--|--------------------|----------------------|-----------------------|---|-------------------|-----------|
| Hotplate stirrer | Ceramic | 7 × 7" | 120 V, 10 A, 1200 W | 310 × 220 × 120 mm (12.2 × 8.8 × 4.8") | 3.5 kg (7.8 lb.) | 76447-040 |
| Hotplate stirrer with NIST traceable calibration certificate | Ceramic | 7 × 7" | 120 V, 10 A, 1200 W | 310 × 220 × 120 mm (12.2 × 8.8 × 4.8") | 3.5 kg (7.8 lb.) | 76447-048 |
| Hotplate stirrer | Ceramic | 10 × 10" | 120 V, 11.2 A, 1344 W | 410 × 280 × 120 mm (16.3 × 11.3 × 4.8") | 5.7 kg (12.6 lb.) | 76447-026 |
| Hotplate stirrer with NIST traceable calibration certificate | Ceramic | 10 × 10" | 120 V, 11.2 A, 1344 W | 410 × 280 × 120 mm (16.3 × 11.3 × 4.8") | 5.7 kg (12.6 lb.) | 76447-028 |



**\$283
VALUE**

Buy a VWR® Professional hotplate stirrer and receive a FREE support rod and clamp kit!

FREE STIRRER SUPPORT ROD AND CLAMP KIT

Use with temperature probes, thermometers or other glass tubes. Includes 18" (45.7 cm) stainless steel support rod, thermometer/temperature probe extension clamp, 3-prong medium swivel clamp and hook connector.

| Length (cm) | Cat. No. |
|-------------|-----------|
| 45.7 | 11301-110 |

Please follow the instructions below to redeem offer.

1. Place your order on vwr.com or call 1-800-932-5000 before December 31, 2025.
2. To redeem offer: have your proof of purchase (VWR invoice or packing slip), visit vwr.com/promotions, search for Promo Code 4292 and fill in the online redemption form by January 15, 2026.

Note: Please specify the block quantity you need for each free block to fit the capacity of your unit. Offer valid only for products listed on this flyer.

Allow 8 weeks for delivery of free items.

Offer valid January 1, 2025 – December 31, 2025.



Avantor
ScienceCentral™
avantorsciences.com

Any prices, product, and/or services details are current when published and subject to change. Offers subject to product availability. Void where prohibited by law or company policy. Any claims concerning products manufactured by third parties are the sole claims of the manufacturer and not those of Avantor or its affiliates. All products are subject to the terms and conditions of sale of Avantor or its affiliates, including the applicable limitations of liability contained therein. The information herein should not be considered a warranty of any kind. Please carefully review any product information sheets and safety information for the product. Trademarks are owned by Avantor Inc. or its affiliates, unless otherwise noted. | Visit vwr.com and avantorsciences.com to view our privacy policy and additional disclaimers. © 2025 Avantor, Inc. All rights reserved.